



Serving Healthcare Worldwide Since 1984



Podiatry Catalogue Nail Nippers and Splitters

Reusable, Non-Sterile, Available in Various
Length & Blade Types

Quality in Stainless Steel

Hand instruments for foot care must satisfy high quality requirements. Only high quality materials, first class craftsmanship and well thought out design guarantee safety and comfort in use.

And: an instrument becomes cost effective and financially viable for the practice if it has a long life and its properties for use remain consistent.

Only instruments made of stainless steel fulfill these basic requirements of the guidelines for medical devices. In contrast to nickel coated instruments, they are entitled to carry the **CE** symbol and are approved as medical devices. With these tested and certified instruments, medical foot care specialists are justified in claiming full compliance with their duty of care regarding hygiene.



Tested for safety

In order to meet the high requirements and expectations of the user in all areas, the **ACUMED** quality requirements are also exceptionally high.

From the onset the company has therefore devised its own quality standard, which stringently restricts the permitted **DIN/ISO** standard tolerances in every respect. This is already noticeable in the selection and testing of the materials and continues throughout the entire production process.

The instrument manufacturing technology based on the symbiosis of high tech and craftsmanship as well as vertical range of production provides compliance with standards and the greatest precision with consistent quality.

At **ACUMED** they have known for more than 40 years that there is nothing more economical than high quality.

The memory
of poor quality
lasts longer
than the brief
pleasure of the
low price.

Recognised Quality

That is how quality is recognised – after long-term use:

- The handles lie comfortably in the hand.
- Scissors and clippers glide smoothly and evenly up and down as the handle is moved.
- Joints are finished so that one half of an instrument fits into the other as it crosses and is then forged again. There is no unnecessary play of movement in the closing area.
- The cutting edges of forceps meet precisely right up to the tip.
- There are no gaps in the scissors between the sharp cutting edges.
- All surfaces are smooth and free of rust or discolouration – even joints, cutting and rubbing surfaces.



**TITANIUM
BLACK**

*Interested in Black Coated
Nail Nippers? Just add (BL)
after the code to order
Titanium Black Coated*

 **ACUMED**
International Co.

Serving Healthcare Worldwide Since 1984





Features & Benefits

Beautifully Packed in Box

Anatomy of Acumed Nail Instruments



JAW DESIGN

Exact blade alignment and precision sharpening facilitates effortless cutting and extends the performance life of the instrument.

BOX LOCK DESIGN

Fully tempered box lock maintains precise instrument alignment.

ELECTRO-CHEMICAL ETCHING

Instruments are etched with the Medcave logo & catalog number for easy re-ordering. This etching technique will not damage the instruments.

SPRINGS

Springs enable controlled jaw movement & reduce hand fatigue associated with some procedures.

PASSIVATION PROCESS

Passivation protects the instrument against corrosion due to repeated sterilization.

ERGONOMIC DESIGN

Outstanding instrument design combines shape, length, weight, balance and feel for optimal comfort and performance.

MIRROR & SATIN FINISH

This finish simplifies the care and cleaning process by resisting the build up of debris thereby reducing corrosion.

STRENGTH & HARDNESS

Manufactured to ensure optimal Rockwell hardness for instrument durability and longer-lasting blade edge sharpness.

Nail Nippers

ACUMED the name of excellence of quality
Over 20 years of experience, we are delivering
top quality chiropody nail nippers worldwide.

Barrel Spring Nail Nipper

Protect against excessivestress on the podiatrist's
hand and elbow joints.



Acumed Reference No.

AP-04-00

Barrel Spring Nail Nipper
15cm



Acumed Reference No.

AP-04-01

Aesculap Reference No.

G495

Barrel Spring Nail Nipper
12cm



Acumed Reference No.

AP-04-02

Aesculap Reference No.

HF250

Barrel Spring Nail Nipper
15cm



Nail Nippers

Nail Nipper

Cutting for the free nail edges.
Strong Model, Particularly Stable.

Acumed Reference No.

AP-04-03

Aesculap Reference No.

HF241

Nail Nipper

11 cm



Acumed Reference No.

AP-04-04

Aesculap Reference No.

HF243

Nail Nipper

13 cm



Acumed Reference No.

AP-04-05

Aesculap Reference No.

HF210

Box Jointed,
Cantilevered Nail Nipper

13.5 cm



Acumed Reference No.

AP-04-06

Aesculap Reference No.

HF211

Box Jointed,
Cantilevered Nail Nipper

13.5 cm



Nail Nipper

Cutting for the free nail edges.

Strong Model, Particularly Stable.



Acumed Reference No.

AP-04-07

Aesculap Reference No.

HF223

Box Joint

Nail Nipper Straight

13 cm



Acumed Reference No.

AP-04-08

Aesculap Reference No.

HF203

Box Joint

Nail Nipper Curved

13 cm

Corner Clipper

Cutting sides and the nail's free distal edges also for the removal of the ingrown area of the nail (Unguis incarnatus) in the nail fold or for the removal of hyperkeratosis and clavi (corns and hard skin).



Acumed Reference No.

AP-04-09

Aesculap Reference No.

HF475

Box Joint

Nail Nipper Straight

11.5 cm



Acumed Reference No.

AP-04-10

Aesculap Reference No.

HF476

Box Joint

Extra Slim Nail Nipper Straight

11.5 cm

Nail Nippers

Corner Clipper

For cutting the sides and the free distal edges of the nails. For removal of the ingrown part of the nail (Unguis incarnatus) in the nail fold or for the removal of hyperkeratosis and clavi (corns and hard skin).

Acumed Reference No.

AP-04-11

Aesculap Reference No.

HF477

Nail Nipper

11.5 cm



Acumed Reference No.

AP-04-12

Aesculap Reference No.

HF481

Nail Nipper

13.5 cm



Acumed Reference No.

AP-04-13

Aesculap Reference No.

HF478

Nail Nipper

13 cm



Acumed Reference No.

AP-04-14

Aesculap Reference No.

HF479

Nail Nipper

13 cm



Acumed Reference No.

AP-04-15

Aesculap Reference No.

HF480

Nail Nipper

15 cm



Corner Clipper

For cutting the sides and the free distal edges of the nails. For removal of the ingrown part of the nail (Unguis incarnatus) in the nail fold or for the removal of hyperkeratosis and clavi (corns and hard skin).

Nail Nippers



Acumed Reference No.

AP-04-16

Aesculap Reference No.

HF483

Nail Nipper

13 cm



Acumed Reference No.

AP-04-17

Aesculap Reference No.

HF482

Nail Nipper

12 cm



Acumed Reference No.

AP-04-18

Aesculap Reference No.

HF490

Nail Nipper

11.5 cm



Acumed Reference No.

AP-04-19

Aesculap Reference No.

HF491

Nail Nipper

10 cm



Acumed Reference No.

AP-04-20

Aesculap Reference No.

HF470

Nail Nipper

10 cm

Nail Nippers

Skin Clipper

For the removal of all kinds of clavi from the feet, such as clavi between the toes, hyperkeratosis and protruding edges of corns.

Acumed Reference No.

AP-04-21

Aesculap Reference No.

HF460

Cantilevered Nail Nippers
10 cm



Acumed Reference No.

AP-04-22

Aesculap Reference No.

HF461

Cantilevered Nail Nippers
11.5 cm



Acumed Reference No.

AP-04-23

Aesculap Reference No.

HF465

Cantilevered Nail Nippers
11.5 cm



Acumed Reference No.

AP-04-24

Aesculap Reference No.

HF448

Nail Nipper
10 cm



Acumed Reference No.

AP-04-25

Aesculap Reference No.

HF466

Nail Splitter
13.5 cm



Nail Nippers

General Purpose Nail Nippers

These Nail Nippers are for general purpose, made with high quality, surgical grade - non rusty stainless steel.



AP-04-26



AP-04-27



AP-04-28



AP-04-29

Reference No.	Size	Description
AP-04-26	14cm	Nail Nipper General Curved B/J Catch
AP-04-27	14cm	Nail Nipper General Straight B/J Catch
AP-04-28	14cm	Nail Nipper General Curved S/J Catch
AP-04-29	14cm	Nail Nipper General Straight S/J Catch



AP-04-30



AP-04-31



AP-04-32



AP-04-33

Reference No.	Acumed	Aesculap	Size	Description
AP-04-30		HF441	8.5cm	Nail Nipper Curved
AP-04-31		HF443	10.5cm	Nail Nipper Curved
AP-04-32		HF445	10.5cm	Nail Nipper Bevelled
AP-04-33		HF446	10.5cm	Nail Nipper Straight

Scissors & Tweezers

Scissors



AP-04-34

Reference No.



AP-04-35



AP-04-36



AP-04-37

Acumed	Aesculap	Size	Description
AP-04-34	HF30	9.5cm	Nail Scissor Curved
AP-04-35	HF31	10.5cm	Nail Scissor Curved
AP-04-36	HF409	9cm	Nail Skin Scissor
AP-04-37	BC324	14.5cm	Dressing Scissor Straight

Tweezers



AP-04-38



AP-04-39



AP-04-40



AP-04-41



AP-04-42



AP-04-43



AP-04-44

Acumed	Aesculap	Size	Description
AP-04-38	OC069	12cm	Adson Anatomical Tweezer
AP-04-39	BD223	15cm	Adson Anatomical Tweezer
AP-04-40	BD303	11.5cm	Splinter Tweezer
AP-04-41	OC220	12cm	Adson Anatomical Tweezer
AP-04-42	BD023	11.5cm	Anatomical Tweezer Wide
AP-04-43	BD964	9.5cm	Splinter Tweezer
AP-04-44	BD312	10.5cm	Splinter Tweezer

Forceps



BEST SELLER

Halstead Mosquito
AP-04-45 - Straight
AP-04-46 - Curved



BEST SELLER

Spencer Wells
AP-04-47 - Straight
AP-04-48 - Curved



BEST SELLER

Blade Removing Forceps
AP-04-49

Reference No.	Size	Description
AP-04-45	12.5cm	Halstead Mosquito Straight Artery Forceps
AP-04-46	12.5cm	Halstead Mosquito Curved Artery Forceps
AP-04-47	14cm	Spencer Wells Straight Artery Forceps
AP-04-48	14cm	Spencer Wells Curved Artery Forceps
AP-04-49	15cm	Blade Removing Forceps



AP-04-50



AP-04-51



AP-04-52



AP-04-53



AP-04-54



AP-04-55



AP-04-56

Reference No.	Size	Description
AP-04-50	0	Gauze Applicator
AP-04-51	128Hz	Tuning Fork
AP-04-52	20cm	Taylor Hammer
AP-04-53	20cm	Dejerine Hammer
AP-04-54	# 3	Scalpel Handle
AP-04-55	# 4	Scalpel Handle
AP-04-56	Standard	Adjustable Scalpel Handle

Elevators and Files



A04-57 **A04-58** **A04-59** **A04-60** **A04-61** **A04-62** **A04-63** **A04-64**
Reference No. **Size** **Description**

A04-57	16.5cm	Mackay Nail Elevator
A04-58	14cm	Fickling Nail Elevator
A04-59	13cm	Fickling Nail Elevator
A04-60	13cm	Fickling Nail Elevator
A04-61	20cm	Diamond Deb
A04-62	1mm	Corner Chisel
A04-63	2mm	Corner Chisel
A04-64	19cm	Lewis Nail Rasp

Nai Probes & Nail Files



Reference # AP-04-65 **Description: Nail Probe Swan Neck - Flat**



Reference # AP-04-66 **Description: Nail Probe Swan Neck - Flat & Fine**



Reference # AP-04-67 **Description: Black Nail File**



Reference # AP-04-68 **Description: Nail Excavator Angled Round Head**



Reference # AP-04-69 **Description: Nail Probe Swan Neck Cup**

Nail Files

Extra Fine Cut



Reference # AP-04-70 Description: Baka Nail File Extra Small Head - Extra Fine Cut

Fine Cut



Reference # AP-04-71 Description: Nail File Medium Head - Fine Cut

Fine Cut



Coarse Cut



Reference # AP-04-72 Description: Double Ended Nail File Medium Heads Fine/Coarse Cut

Curved Head



Straight Head



Reference # AP-04-73 Description: Double Ended Nail File Narrow Head

Nail File



Excavator



Reference # AP-04-74 Description: Nail Excavator & File - Double Ended

Nail File



Reference # AP-04-75 Description: Double Ended Nail File Extra Fine - Medium Heads

Nail File



Reference # AP-04-76 Description: Nail File Swan Neck Probe - Double Ended

Straight Head



Angled Head



Reference # AP-04-77 Description: Double Ended Nail File Medium Heads - Medium Cut

Straight Head



Reference # AP-04-78 Description: Nail File Small Head - Medium Cut



Files, Elevators, Probes,

Medium Cut



Reference #AP-04-79 Description: Nail File Medium Head - Medium Cut - Stainless Steel Made

Fine Cut



Reference #AP-04-80 Description: Nail File Extra Narrow Head - Fine Cut - Stainless Steel Made

1mm Scoop with hole



2mm Scoop with hole



Reference #AP-04-81 Description: Double Ended Nail Curette 1mm & 2mm Scoop with hole - Stainless Steel

2mm Scoop with hole



3mm Scoop with hole



Reference #AP-04-82 Description: Double Ended Nail Curette 2mm & 3mm Scoop with hole - Stainless Steel



Reference #AP-04-83 Description: MacDonald Dissector and Elevator- Double Ended - Stainless Steel Made

Round



Straight



Reference #AP-04-84 Description: Swan Neck Nail Elevator - Double Ended - Stainless Steel Made



Reference #AP-04-85 Description: Nail Spatula - Double Ended - Stainless Steel Made



Reference #AP-04-86 Description: Nail Spatula and Probe - Stainless Steel

Angled



Reference #AP-04-87 Description: Nail Probe Flat Head - Stainless Steel Made



Acumed Logo and Name are Registered and Trade Mark of Acumed



Scan our mobile website to your smartphone for quick reference

Acumed International Co.
Corporate Head
Jinnah Town, Capital Road,
Sialkot - 51310, Pakistan.

Email: sales@acumedco.com
Website: www.acumedco.com