

# Dental Instruments



# AD

Diagnostics

Diagnóstico

Diagnose

*Quality makes Dentistry Better*

 **ACUMED**  
International Co.

Serving Healthcare Worldwide Since 1984

## Know your stems

### Simple Stem

Threading diameter is same as shaft

Consistant diameter shaft

European Pattern



### Cone Socket

Threading diameter is smaller than shaft

Beveled flae locks head into handle

American Pattern



It is important to match a simple stem mirror head to a simple stem mirror handle & cone socket mirror head to cone socket mirror handle



The stream-lined and ergonomic design of the double-sided mirror allows for maximum patient conform and efficient cleaning due to crevice-free construction.



European



American



### Simple Stem

- ADM07S
- #4DS
- Ø 22mm
- ADM08S
- #5DS
- Ø 24mm

### Cone Socket

- ADM07C
- #4DS
- Ø 22mm
- ADM08C
- #5DS
- Ø 24mm



### Simple Stem

- ADM09
- Ø 22mm
- ADM10
- Ø 24mm

### Cone Socket

- ADM11
- Ø 22mm
- ADM12
- Ø 24mm

Simple Stem Pack of 12 Plane	Cone Socket Pack of 12 Plane
Code # 3 - ADM01P	Code # 3 - ADM04P
Code # 4 - ADM02P	Code # 4 - ADM05P
Code # 5 - ADM03P	Code # 5 - ADM06P
Magnifying	Magnifying
Code # 3 - ADM01M	Code # 3 - ADM04M
Code # 4 - ADM02M	Code # 4 - ADM05M
Code # 5 - ADM03M	Code # 5 - ADM06M
Front Surface	Front Surface
Code # 3 - ADM01F	Code # 3 - ADM04F
Code # 4 - ADM02F	Code # 4 - ADM05F
Code # 5 - ADM03F	Code # 5 - ADM06F



Laryngeal Mirrors	
AD-LM1	Ø 10mm
AD-LM2	Ø 12mm
AD-LM3	Ø 14mm
AD-LM4	Ø 16mm
AD-LM5	Ø 18mm
AD-LM6	Ø 20mm
AD-LM7	Ø 22mm
AD-LM8	Ø 24mm
AD-LM9	Ø 26mm
AD-LM10	Ø 28mm
AD-LM11	Ø 30mm

### Laryngeal Mirrors Handle



### Simple Stem

- ADM13
- #4DS
- Ø 22mm
- ADM14
- #5DS
- Ø 24mm



### Cone Socket

- ADM15
- #4DS
- Ø 22mm
- ADM16
- #5DS
- Ø 24mm



### European Front Surface Mouth Mirrors

Double sided simple stem mirror used for indirect vision, retraction of lips, cheks, and tongue, and illumination of dark areas. Available in both European (Simple Stem) American Thread (Cone Socket) threads.



<b>ADMD4</b> Double Sided - Rhodium N.4 Ø 22mm	<b>ADMD5</b> Double Sided - Rhodium N.5 Ø 24mm	<b>ADMD4B</b> Double Sided - Rhodium N.4 Ø 22mm Titanium Black Coating	<b>ADMD5B</b> Double Sided - Rhodium N.5 Ø 24mm Titanium Black Coating
--	--	---	---

## Premium Handle

### Feather Light

Available in Different Colours



## Elite Handle

ULTRA-LIGHT  
Stainless Steel Made

Available in Different Colours



**NEW**

Medcave's newest and most advanced technically designed handle! 9.5mm diameter as our popular round designed handle. More gradual taper from the tip for more comfortable grip even more closer to the tip. Silicon Color grip gives you maximum slip-free gripping power. The handle is desinged with high grade stainless steel to make it more easy clear and sterilize by any method. 14 color rings are available for your choice...



AD-01 (Cone Socket) AD-02 (Simple Stem)



AD-03 (Cone Socket) AD-04 (Simple Stem)



AD-05 (Cone Socket) AD-06 (Simple Stem)



**NEW**

AD-07 (Cone Socket) AD-08 (Simple Stem)

The Perfect Comfort round handle, New developed innovative design handle by our expert engineers to meet the Dentist's needs. 9mm diameter hollow handle is designed for you with a very thin wall, which is very light weight, in order to enable you more comfortable, more ergonomic and more impact working. Its special grip profile with silicon grip and sand blast surface grants the handle Unique Look. Sand Blasted surface minimizes annoying light reflections during working with loupe and microscope.



AD-09 (Cone Socket) AD-10 (Simple Stem)



AD-11 (Cone Socket) AD-12 (Simple Stem)



AD-13 (Cone Socket) AD-14 (Simple Stem)



AD-15 (Cone Socket) AD-16 (Simple Stem)



AD-17 (Cone Socket) AD-18 (Simple Stem)

### Extra Comfort Silicon Grip Handles

**COMFORT GRIP** Set of 5 Handles  
AD-CG5



AD-19 (Cone Socket) AD-20 (Simple Stem)



AD-21 (Cone Socket) AD-22 (Simple Stem)



AD-23 (Cone Socket) AD-24 (Simple Stem)



AD-25 (Cone Socket) AD-26 (Simple Stem)



AD-27 (Cone Socket) AD-28 (Simple Stem)

ALUMINUM Handles Set of 5 Handles  
AG-AH5



AD-29 (Cone Socket) AD-30 (Simple Stem)




AD-31 (Cone Socket) AD-32 (Simple Stem)



AD-33 (Cone Socket) AD-34 (Simple Stem)



AD-35 (Cone Socket) AD-36 (Simple Stem)



AD-37 (Cone Socket) AD-38 (Simple Stem)

All Acumed Handles are made of High Grade Stainless Steel Handles can be autoclaved, chemiclaved, ultrasonically cleaned & sterlized by cold or dry heat method

Cone Socket - AD-39



Simple Stem - AD-40

Small Round Mirror Handle, Cone Socket 12 cm

Cone Socket - AD-41



Simple Stem - AD-42

Small Round Mirror Handle, Cone Socket, 60 mm - 12 cm

Cone Socket - AD-43



Simple Stem - AD-44

Standard Round Handle, 6mm diameter with 30 grams weight

Cone Socket - AD45



Simple Stem - AD-46

Small Octagonal Mirror Handle, 12.5 cm

Cone Socket - AD-47



Hand Form Mirror Handle

Cone Socket - AD-49



Simple Stem - AD-50

Commonly Used Mirror Handle

Cone Socket - AD-51



Simple Stem - AD-52

Cone Socket - AD-53



Simple Stem - AD-54

Commonly Used Mirror Handle

Cone Socket - AD-55



Simple Stem - AD-56

Sand Blast Finish 8mm Hollow Handle



### Handle Advantages:

- Light Weight and very ergonomic handling
- No Problem for patients with allergic reaction
- Attractive Color for childrens and women



Set of 5 Pcs  
AD-SGHS

Simple Stem - AD-58



Cone Socket - AD-57

Simple Stem - AD-58TBK



Cone Socket - AD-57TBK

Simple Stem - AD-60



Cone Socket - AD-59

Acumed Specially Designed Endo Mirror Handles, with Endodontic Scale

Simple Stem - AD-60TBK



Cone Socket - AD-59TBK

Simple Stem - AD-62



Cone Socket - AD-61

Simple Stem - AD-64



Cone Socket - AD-63

Simple Stem - AD-66



Cone Socket - AD-65

Simple Stem - AD-68



Cone Socket - AD-67

Simple Stem - AD-70



Cone Socket - AD-69

**COTTON AND DRESSING TWEEZERS**



Serrated

Grooved,  
Suited to hold silver pins,  
nerve-broaches and for  
introducing medicines.

<b>Dressing</b> 15cm AD-71	<b>Meriam</b> 16cm AD-72	<b>College</b> 15cm AD-73	<b>Perry</b> 13cm AD-74	<b>Castroviejo</b> Lock 15cm AD-75	<b>Castroviejo</b> Lock 15cm AD-76
----------------------------------	--------------------------------	---------------------------------	-------------------------------	---	---

**SERRATED POINTS - LIGHT WEIGHT FOR BETTER HANDLING**

**SUTURE TWEEZERS**



<b>Perry</b> 13cm AD-77	<b>Flag</b> 15cm AD-78	<b>College</b> Lock 15cm AD-79	<b>Meriam</b> 16cm AD-80	<b>College</b> 15cm AD-81	<b>College</b> 15cm AD-82	<b>Suture</b> 15.5cm AD-83	<b>Suture</b> 15.5cm AD-84
-------------------------------	------------------------------	---	--------------------------------	---------------------------------	---------------------------------	----------------------------------	----------------------------------

All our explorers are also available in 6mm solid handle and 8mm Titanium Black handles.

## Classic

### DID YOU KNOW?

The use of explorers is to examine teeth for decay (carries), calculus and other abnormalities. Explorers have flexible thin tips that taper to a sharp point and come in various single end and double end patterns. All Acumed Explorer tips are made from spring tempered stainless steel and are machine ground exacting tolerances. Each tip is hand sharpened and formed into its final shape by expert craftsmen. With Acumed, you can assured of the premium quality, consistency, and the sharpest tips in the industry.

*Quality makes Dentistry Better*



AD-85 Fig. 2	AD-86 Fig. 23	AD-87 Fig. 9	AD-88 Fig. 16	AD-89 Fig. 11	AD-90 Fig. 3	AD-91 Fig. 3A	AD-92 Fig. 5
-----------------	------------------	-----------------	------------------	------------------	-----------------	------------------	-----------------



AD-93 Fig. 6	AD-94 Fig. 6A	AD-95 Fig. 8A	AD-96 Fig. 10	AD-97 Fig. 12	AD-98 Fig. 17	AD-99 Fig. 17A	AD-100 Fig. 18
-----------------	------------------	------------------	------------------	------------------	------------------	-------------------	-------------------


*Matt Finish  
8mm Hollow Handle  
Fully Autoclavable*

Single Ended Dental Explorers


# Classic

*Matt Finish  
8mm Hollow Handle  
'Fully Autoclavable'*

# Diagnostics



Periodontal probes with Titanium Black Coating and Laser Marked Black background to highlight the contrast ratio for better reading  
Rounded tip to prevent traumatic issue  
Accurate millimeter scale markings realized by laser  
Atraumatic scale non-grooved for a smooth transition




Single Ended Dental Probes

## Classic


*Matt Finish  
8mm Hollow Handle  
Fully Autoclavable*



**AD-122**  
Fig. 23/17



**AD-123**  
Fig. 23/17A



**AD-124**  
Fig. 23/6



**AD-125**  
Fig. 2



**AD-126**  
Fig. 3



**AD-127**  
Fig. 3/6



**AD-128**  
Fig. 5/17A



**AD-129**  
Fig. 9/17



**AD-130**  
Fig. 11/12

### DID YOU KNOW?

ACUMED has top quality Dental exploration kits, please go to our page number 11 of diagnostic section. ACUMED has different types of dental exploration kit to provide best quality kits for doctors, for better examen during diagnostics.

*Quality makes Dentistry Better*





# Classic

*Matt Finish  
8mm Hollow Handle  
Fully Autoclavable*

# Diagnostics



AD-131 Fig. 3/Williams



AD-132 Fig. 3/UM



AD-133 Fig. 3/Goldman Fox



AD-134 Fig. 16/Goldman Fox



AD-135 Fig. 16/Williams



AD-136 Fig. 16/UM



AD-137 Fig. 17/Goldman Fox



AD-138 Fig. 17/Williams



AD-139 Fig. 17/UM




AD-140 Fig. 17/WHO



AD-141 Fig. 17/CP12



AD-142 Fig. 23/WHO



AD-143 Fig. 23/CP12



AD-144 Fig. 23/UM



AD-145 Fig. 23/Williams



AD-146 Fig. 23/Goldman Fox

6mm Solid Handles



**AD-147**

- 1 Pc - Merriam Tweezer
- 1 Pc - Mouth Mirror # 4
- 1 Pc - Mirror Handle
- 1 Pc - Single Ended Probe # 9

**AD-148**

- 1 Pc - Merriam Tweezer
- 1 Pc - Mouth Mirror # 4
- 1 Pc - Mirror Handle
- 1 Pc - Single Ended Probe # 23

**AD-149**

- 1 Pc - College Tweezer
- 1 Pc - Mouth Mirror # 4
- 1 Pc - Mirror Handle
- 1 Pc - Double Ended Probe # 5

**AD-150**

- 1 Pc - Merriam Tweezer
- 1 Pc - Mouth Mirror # 4
- 1 Pc - Mirror Handle
- 1 Pc - Single Ended Probe # 23



*Classic*

ALUMINUM Handles



**AD-151**

- 1 Pc - Merriam Tweezer
- 1 Pc - Mouth Mirror # 4
- 1 Pc - Silicon Mirror Handle
- 1 Pc - Silicon Single Ended Probe # 23
- 1 Pc - Silicon CP 15 Probe
- 1 Pc - Plastic Sterilization Tray



**AD-152**

- 1 Pc - Merriam Tweezer
- 1 Pc - Mouth Mirror # 4
- 1 Pc - 8mm Hollow Mirror Handle
- 1 Pc - 8mm Handle Single Ended Probe # 23
- 1 Pc - 8mm Handle Goldman Fox Probe
- 1 Pc - Stainless Steel Sterilization Tray



**AD-153**

- 1 Pc - Merriam Tweezer
- 1 Pc - Mouth Mirror # 4
- 1 Pc - Mirror Handle
- 1 Pc - Single Ended Probe # 23
- 1 Pc - WHO Probe

## Dental Photography Mirrors



<b>AD-154</b> Occlusal (Child)	<b>AD-155</b> Occlusal (Child Long)	<b>AD-156</b> Occlusal (Adult)	<b>AD-157</b> Occlusal (Adult Long)	<b>AD-158</b> Occlusal (Lingual)	<b>AD-159</b> Occlusal (Buccal Narrow)	<b>AD-160</b> Occlusal (Buccal)
--------------------------------------	---	--------------------------------------	---	--	--	---------------------------------------



<b>AD-161</b> Child Occlusal Adult Occlusal	<b>AD-162</b> Adult Buccal Adult Occlusal	<b>AD-163</b> Child Occlusal Adult Occlusal	<b>AD-164</b> Child Occlusal Child Buccal	<b>AD-165</b> Lingual Upper Lingual Lower	<b>AD-166</b> Adult Buccal Child Buccal	<b>AD-167</b> Buccal Buccal
---	---	---	---	---	---	-----------------------------------

## Dental Photography Image Contrast Retractors Black - Set of 6 Pcs

Code: AD-PGRS





*Acumed Logo and Name are Registered and Trade Mark of Acumed*



*Scan our mobile website to your smartphone for quick reference*

**Acumed International Co.**  
**Corporate Head**  
**Jinnah Town, Capital Road,**  
**Sialkot - 51310, Pakistan.**

Email: [sales@acumedco.com](mailto:sales@acumedco.com)

Website: [www.acumedco.com](http://www.acumedco.com)